## **CAPSULAR TENSION RINGS**

## Flexible Insertion, Stable Result

Highly flexible due to the compression molding technology. These PMMA Capsular Tension Rings are indicated for the stabilization of weakened, broken or missing zonules preventing IOL decentration after capsular shrinking.



## Capsular Tension Ring | Model 275 12/10



Compression: From 12 mm to 10 mm

Material: PMMA | Flexible Ring

Available in: Clear, brown and green PMMA

Overall Ø: 12 mm

## Capsular Tension Ring | Model 276 13/11



Compression: From 13 mm to 11 mm

Material: PMMA | Flexible Ring

Available in: Clear, brown and green PMMA

Overall Ø: 13 mm

- ✓ CIRCULAR EXPANSION AND

  STABILIZATION OF THE CAPSULAR BAG
- ✓ SAFE IOL CENTRATION IN EYES WITH ZONULAR DEHISCENCE
- ✓ PREVENTS IOL DECENTRATION AFTER CAPSULAR SHRINKING
- ✓ STABILIZED CONDITIONS DURING
  (PHACO EMULSIFICATION) SURGERY
- **✓ REDUCED RISK OF CAPSULAR FIBROSIS**



Recommended to be used in combination with the ON 385 EASY Control™ Micro Inserter

www.ophtec.com